



DUO-FLOW®


Indications for Use: The Medcomp® Duo-Flow® Catheter is indicated for use in attaining Short-Term vascular access for Hemodialysis and Apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include subclavian vein or femoral vein as required. The curved Duo-Flow® Catheter is intended for internal jugular vein insertion. This catheter is indicated for a duration less than (30) days. For femoral placement, monitor catheter condition closely.

Contraindications: This catheter is intended for Short-Term vascular access only and should not be used for any purpose other than indicated in these instructions.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



Polyurethane Material	Tapered Tip Tip Design	9F, 11.5F French Sizes	Straight, Pre-Curved, Raulerson IJ Configurations
---------------------------------	----------------------------------	----------------------------------	---



TAPERED TIP



COAXIAL INTERNAL LUMEN DESIGN



RAULERSON IJ



CATHETER ONLY		
XTP114CT	11.5F X 12CM DUO-FLOW®, STRAIGHT	5/BOX
XTP116CT	11.5F X 15CM DUO-FLOW®, STRAIGHT	5/BOX
XTP118CT	11.5F X 20CM DUO-FLOW®, STRAIGHT	5/BOX
XTP119CT	11.5F X 24CM DUO-FLOW®, STRAIGHT	5/BOX
XTP94IJC	9F X 12CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP96IJC	9F X 15CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP98IJC	9F X 20CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP114IJC	11.5F X 12CM, DUO-FLOW®, RAULERSON IJ	5/BOX
XTP116IJC	11.5F X 15CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP118IJC	11.5F X 20CM DUO-FLOW®, RAULERSON IJ	5/BOX

CATHETER ONLY SET
CONTENTS:
(1) Catheter
(1) Dilator
(2) End Caps

BASIC SET		
XTP94MT	9F X 12CM DUO-FLOW®, STRAIGHT	5/BOX
XTP96MT	9F X 15CM DUO-FLOW®, STRAIGHT	5/BOX
XTP98MT	9F X 20CM DUO-FLOW®, STRAIGHT	5/BOX
XTP114MTE	11.5F X 12CM DUO-FLOW®, STRAIGHT	5/BOX
XTP116MTE	11.5F X 15CM DUO-FLOW®, STRAIGHT	5/BOX
XTP118MTE	11.5F X 20CM DUO-FLOW®, STRAIGHT	5/BOX
XTP119MTE	11.5F X 24CM DUO-FLOW®, STRAIGHT	5/BOX
XTP114IJS-2	11.5F X 12CM DUO-FLOW®, PRE-CURVED	5/BOX
XTP116IJS-2	11.5F X 15CM DUO-FLOW®, PRE-CURVED	5/BOX
XTP118IJS-2	11.5F X 20CM DUO-FLOW®, PRE-CURVED	5/BOX
XTP94IJS	9F X 12CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP96IJS	9F X 15CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP98IJS	9F X 20CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP114IJS-2	11.5F X 12CM, DUO-FLOW®, RAULERSON IJ	5/BOX
XTP115IJS-2	11.5F X 13.5CM, DUO-FLOW®, RAULERSON IJ	5/BOX
XTP116IJS-2	11.5F X 15CM DUO-FLOW®, RAULERSON IJ	5/BOX
XTP118IJS-2	11.5F X 20CM DUO-FLOW®, RAULERSON IJ	5/BOX

BASIC SET CONTENTS:
(1) Catheter
(1) Introducer Needle
(1) Guidewire
(1) Dilator
(1) Scalpel*
(1) Adhesive Wound Dressing
(2) End Caps**

*9F Raulerson IJ Only
**9F and 11.5F Pre-Curved Basic Set Only

Refer to the Table of Contents for Important Risk Information regarding this device.